







SAMAY
Surgical Pvt. Ltd.

The Ultimate Orthopedic Implant Solution



MAXILLOFACIAL PLATES
SYSTEM

www.samaysurgical.com

Introduction

Vision

Our ultimate vision is to strong provider of orthopedic implants for surgeons and hospitals around the world, known for our commitment to excellence and customer satisfaction.



Mission

To providing safe, effective and affordable orthopedic implant with Marvelous quality.



Our organization, Samay Surgical, derives its name from the Sanskrit word 'samay' (समय), which means the appointed or proper time - the right moment for doing anything. This philosophy is reflected in our work, as we believe in investing time up-front in developing systems and processes that promote long-term growth and sustainability.


We understand that the needs and expectations of our customers vary depending on their location, culture, and healthcare system. That's why we take a customer-centric approach to our operations, focusing on understanding and meeting the unique needs of each market we serve.

Celebrating over a decade of growth and progress, Samay Surgical has come a long way since its inception in 2006 as a small company that only produced screws. Today, our product range has expanded significantly, and we are proud to offer a diverse portfolio that includes plates, nails, wires, pins, spine, foot systems, ankle systems, mini-micro plate systems, and much more.

However, our commitment to excellence does not stop with just our products. We are also dedicated to providing the best possible customer service, enhancing the customer experience, and cultivating a positive company culture. At Samay Surgical, we believe that it is possible to simultaneously deliver the ultimate solution and cost-effective goods to our customers, inspiring the world with our success.

Through consistent organic growth and strategic acquisitions, we have transformed ourselves from a small manufacturing plant into the largest orthopedic implant manufacturer in the industry, supporting the orthopedic needs of people all around the world.



At Samay Surgical, we are passionate about our mission to provide the highest quality products, exceptional customer service, and an unparalleled customer experience. We are committed to continuing our growth and innovation, always striving to be at the forefront of our industry and inspiring others to achieve excellence in their own endeavors.

Mini Plate Straight 1.5mm (Without Gap)

TT Gr.5	Size
T2197.02	2 Hole
T2197.03	3 Hole
T2197.04	4 Hole
T2197.05	5 Hole
T2197.06	6 Hole
T2197.08	8 Hole
T2197.10	10 Hole
T2197.12	12 Hole
T2197.16	16 Hole
T2197.20	20 Hole
T2197.30	30 Hole
T2197.50	50 Hole

Mini Plate Straight 1.5mm (With Gap)

TT Gr.5	Size
T2198.02	2 Hole
T2198.03	3 Hole
T2198.04	4 Hole
T2198.05	5 Hole
T2198.06	6 Hole
T2198.08	8 Hole



**Mini 'L' Shape 90° Plate 1.5mm(Without Gap)**

Left Right

	TT Gr.5	Size
Right	T2199R.22	2+2 Hole
Left	T2199L.22	2+2 Hole

Mini 'L' Shape 90° Plate Right 1.5mm


TT Gr.5	Size
T2200R.22	2+2 Hole
T2200R.23	2+3 Hole
T2200R.32	3+2 Hole
T2200R.33	3+3 Hole

**Mini 'L' Shape 90° Plate Left 1.5mm**


TT Gr.5	Size
T2200L.22	2+2 Hole
T2200L.23	2+3 Hole
T2200L.32	3+2 Hole
T2200L.33	3+3 Hole

Mini 'L' Shape 100° Plate 1.5mm

	TT Gr.5	Size
Right	T2201R.32	3+2 Hole
Left	T2201L.32	3+2 Hole





Mini 'L' Shape 120° Plate 1.5mm (Without Gap)



	TT Gr.5	Size
Right	T2202R.22	2+2 Hole
Left	T2202L.22	2+2 Hole

Mini 'L' Shape 120° Plate Right 1.5mm

TT Gr.5	Size
T2203R.22	2+2 Hole
T2203R.23	2+3 Hole
T2203R.32	3+2 Hole
T2203R.33	3+3 Hole




Mini 'L' Shape 120° Plate Left 1.5mm

TT Gr.5	Size
T2203L.22	2+2 Hole
T2203L.23	2+3 Hole
T2203L.32	3+2 Hole
T2203L.33	3+3 Hole

Mini Plate Straight 1.5mm (With Long Gap)

TT Gr.5	Size
T2291.02	2 Hole




Mini 'T' Shape Plate 1.5mm


TT Gr.5	Size
T2204.22	2+2 Hole
T2204.32	3+2 Hole
T2204.33	3+3 Hole

Mini 'Y' Shape Plate 1.5mm (Without Gap)


TT Gr.5	Size
T2205.22	2+2 Hole
T2205.23	2+3 Hole

**Mini 'Y' Shape Plate 1.5mm (With Gap)**

TT Gr.5	Size
T2206.22	2+2 Hole
T2206.23	2+3 Hole

**Mini Double 'Y' Shape Plate 1.5mm**

TT Gr.5
T2207.00







Mini 'I' Shape Plate 1.5mm

TT Gr.5
T2208.00

Mini 'H' Shape Plate 1.5mm

TT Gr.5
T2209.00





Mini 'X' Shape Plate 1.5mm

TT Gr.5
T2210.00


Mini Angular Plate 150° Right 1.5mm

TT Gr.5	Size
T2211R.04	4 Hole
T2211R.06	6 Hole



Mini Angular Plate 150° Left 1.5mm


TT Gr.5	Size
T2211L.04	4 Hole
T2211L.06	6 Hole

Mini Curved Plate Right 1.5mm

TT Gr.5	Size
T2212R.04	4 Hole
T2212R.06	6 Hole

Mini Curved Plate Left 1.5mm

TT Gr.5	Size
T2212L.04	4 Hole
T2212L.06	6 Hole

**Mini 'Z' Shape Plate Type-1 1.5mm**

Left




Right


TT Gr.5	Size
T2213R.00	Right
T2213L.00	Left

Mini 'Z' Shape Plate Type-2 1.5mm

TT Gr.5	Size
T2214R.00	Right
T2214L.00	Left




Left



Right


Mini Orbita Plate 1.5mm (With Gap)




TT Gr.5	Size
T2215.04	4 Hole
T2215.06	6 Hole
T2215.08	8 Hole

Mini Orbita Plate 1.5mm (Without Gap)

TT Gr.5	Size
T2216.04	4 Hole
T2216.06	6 Hole
T2216.08	8 Hole
T2216.10	10 Hole
T2216.12	12 Hole
T2216.14	14 Hole
T2216.16	16 Hole
T2216.18	18 Hole
T2216.20	20 Hole




Mini Orbital Floor Plate 1.5mm




TT Gr.5	Size
T2217.00	Small
T2218.00	Medium
T2219.00	Big


Mini Nose Plate 1.5mm

TT Gr.5	Size
T2220.00	H-Shape
T2221.00	Y-Shape




Mini Chin Plate 1.5mm


TT Gr.5	Size
T2222.08	8.0mm
T2222.10	10mm

**Mini Burr Holes Prebend Plate**




TT Gr.5	Size
T2305.18	18mm
T2305.22	22mm

**Mini Burr Holes Drainage Plate**

TT Gr.5	Size
T2296.18	18mm
T2296.22	22mm

**Micro Burr Holes Plate 1.5mm**




TT Gr.5	Size
T2223.12	12mm

Mini 3D Plate 1.5mm


TT Gr.5	Size
T2224.22	2+2 Hole
T2224.23	2+3 Hole
T2224.24	2+4 Hole

Mini Square Shape Plate 1.5mm

TT Gr.5	Size
T2225.22	2+2 Hole
T2225.23	2+3 Hole
T2225.24	2+4 Hole





Eazy Mesh 1.5mm



TT Gr.5	Size
T2226.11	1' x 1'
T2226.22	2' x 2'
T2226.23	2' x 3'
T2226.24	2' x 4'
T2226.26	2' x 6'
T2226.33	3' x 3'
T2226.44	4' x 4'
T2226.46	4' x 6'
T2226.55	5' x 5'
T2226.66	6' x 6'
T2226.86	8' x 6'


Wire Mesh 1.5mm

TT Gr.5	Size
T2228.11	1' x 1'
T2228.22	2' x 2'
T2228.23	2' x 3'
T2228.24	2' x 4'
T2228.33	3' x 3'
T2228.44	4' x 4'
T2228.46	4' x 6'
T2228.48	4' x 8'



Eazy Mesh Pre Bend 1.5mm



TT Gr.5	Size
T2227.11	1' x 1'
T2227.22	2' x 2'
T2227.23	2' x 3'
T2227.24	2' x 4'
T2227.26	2' x 6'
T2227.33	3' x 3'
T2227.44	4' x 4'
T2227.46	4' x 6'
T2227.55	5' x 5'
T2227.66	6' x 6'
T2227.86	8' x 6'

Mini Plate Straight 2.0mm (Without Gap)

TT Gr.5	Size
T2229.02	2 Hole
T2229.03	3 Hole
T2229.04	4 Hole
T2229.05	5 Hole
T2229.06	6 Hole
T2229.08	8 Hole
T2229.10	10 Hole
T2229.12	12 Hole
T2229.16	16 Hole
T2229.20	20 Hole
T2229.30	30 Hole
T2229.50	50 Hole

Mini Plate Straight 2.0mm (With Gap)



TT Gr.5	Size
T2230.02	2 Hole
T2230.03	3 Hole
T2230.04	4 Hole
T2230.05	5 Hole
T2230.06	6 Hole
T2230.08	8 Hole

**Mini 'L' Shape 90° Plate Right 2.0mm(Without Gap)**


	TT Gr.5	Size
Right	T2231R.22	2+2 Hole
Left	T2231L.22	2+2 Hole

Mini 'L' Shape 90° Plate Right 2.0mm

TT Gr.5	Size
T2232R.22	2+2 Hole
T2232R.23	2+3 Hole
T2232R.32	3+2 Hole
T2232R.33	3+3 Hole

**Mini 'L' Shape 90° Plate Left 2.0mm**

TT Gr.5	Size
T2232L.22	2+2 Hole
T2232L.23	2+3 Hole
T2232L.32	3+2 Hole
T2232L.33	3+3 Hole





Mini 'L' Shape 100° Plate 2.0mm

	TT Gr.5	Size
Right	T2233R.32	3+2 Hole
Left	T2233L.32	3+2 Hole

Mini 'L' Shape 120° Plate Right 2.0mm(Without Gap)

	TT Gr.5	Size
Right	T2234R.22	2+2 Hole
Left	T2234L.22	2+2 Hole




Mini 'L' Shape 120° Plate Right 2.0mm

	TT Gr.5	Size
	T2235R.22	2+2 Hole
	T2235R.23	2+3 Hole
	T2235R.32	3+2 Hole
	T2235R.33	3+3 Hole


Mini 'L' Shape 120° Plate Left 2.0mm

TT Gr.5	Size
T2235L.22	2+2 Hole
T2235L.23	2+3 Hole
T2235L.32	3+2 Hole
T2235L.33	3+3 Hole



Mini Plate Straight 2.0mm (With Long Gap)


TT Gr.5	Size
T2292.02	2 Hole

Mini 'T' Shape Plate 2.0mm


TT Gr.5	Size
T2236.22	2+2 Hole
T2236.32	3+2 Hole
T2236.33	3+3 Hole

Mini 'Y' Shape Plate 2.0mm (Without Gap)


TT Gr.5	Size
T2237.22	2+2 Hole
T2237.23	2+3 Hole

**Mini 'Y' Shape Plate 2.0mm (With Gap)**

TT Gr.5	Size
T2238.22	2+2 Hole
T2238.23	2+3 Hole

**Mini Double 'Y' Shape Plate 2.0mm**

TT Gr.5
T2239.00



Mini 'I' Shape Plate 2.0mm

TT Gr.5


T2240.00



Mini 'H' Shape Plate 2.0mm

TT Gr.5


T2241.00



Mini 'X' Shape Plate 2.0mm

TT Gr.5

T2242.00



Mini Angular Plate 150° Right 2.0mm

TT Gr.5

Size

T2243R.04 4 Hole

T2243R.06 6 Hole



Mini Angular Plate 150° Left 2.0mm

TT Gr.5

Size


T2243L.04 4 Hole

T2243L.06 6 Hole




Mini Curved Plate Right 2.0mm


TT Gr.5	Size
T2244R.04	4 Hole
T2244R.06	6 Hole

**Mini Curved Plate Left 2.0mm**

TT Gr.5	Size
T2244L.04	4 Hole
T2244L.06	6 Hole

**Mini 'Z' Shape Plate Type-I 2.0mm**

TT Gr.5	Size
T2245R.00	Right
T2245L.00	Left




Left

Right

Mini 'Z' Shape Plate Type-2 2.0mm


TT Gr.5	Size
T2246R.00	Right
T2246L.00	Left



Left

Right


Mini Orbita Plate 2.0mm (With Gap)



TT Gr.5	Size
T2247.04	4 Hole
T2247.06	6 Hole
T2247.08	8 Hole

Mini Orbita Plate 2.0mm (Without Gap)

TT Gr.5	Size
T2248.04	4 Hole
T2248.06	6 Hole
T2248.08	8 Hole
T2248.10	10 Hole
T2248.12	12 Hole
T2248.14	14 Hole
T2248.16	16 Hole
T2248.18	18 Hole
T2248.20	20 Hole




Mini Chin Plate 4 Hole 2.0mm



TT Gr.5	Size
T2249.08	8mm
T2249.10	10mm


Micro Burr Holes Plate 2.0mm

TT Gr.5	Size
T2250.00	12mm




Mini Burr Holes Prebend Plate

TT Gr.5	Size
T2306.18	18mm
T2306.22	22mm



Mini Burr Holes Drainage Plate

TT Gr.5	Size
T2297.18	18mm
T2297.22	22mm



Mini Nose Plate

TT Gr.5	Size
T2251.00	H-Shape
T2252.00	Y-Shape

Mini Twisted Plate 2.0mm

	TT Gr.5	Size
Right	T2253R.05	5 Hole
Left	T2253L.05	5 Hole




Mini 3D Plate 2.0mm

TT Gr.5	Size
T2254.22	2+2 Hole
T2254.23	2+3 Hole
T2254.24	2+4 Hole

Mini Square Shape Plate 2.0mm

TT Gr.5	Size
T2255.22	2+2 Hole
T2255.23	2+3 Hole
T2255.24	2+4 Hole



Mini Delta Plate 2.0mm

TT Gr.5


T2256.00



Mini Grid Plate 2.0mm

TT Gr.5

T2257.00



Mini Trapezoid Condyle Plate 2.0mm

TT Gr.5

T2258.00



Mini Chin Plate 5 Hole 2.0mm

TT Gr.5


Size


T2259.06 6 Hole

T2259.08 8 Hole

T2259.10 10 Hole

T2259.12 12 Hole





Eazy Mesh 2.0mm

TT Gr.5	Size
T2260.11	1' x 1'
T2260.22	2' x 2'
T2260.23	2' x 3'
T2260.24	2' x 4'
T2260.26	2' x 6'
T2260.33	3' x 3'
T2260.44	4' x 4'
T2260.46	4' x 6'
T2260.55	5' x 5'
T2260.66	6' x 6'
T2260.86	8' x 6'


Eazy Mesh Pre Bend 2.0mm

TT Gr.5	Size
T2261.11	1' x 1'
T2261.22	2' x 2'
T2261.23	2' x 3'
T2261.24	2' x 4'
T2261.26	2' x 6'
T2261.33	3' x 3'
T2261.44	4' x 4'
T2261.46	4' x 6'
T2261.55	5' x 5'
T2261.66	6' x 6'
T2261.86	8' x 6'

**Wire Mesh 2.0mm**

TT Gr.5	Size
T2262.11	1' x 1'
T2262.22	2' x 2'
T2262.23	2' x 3'
T2262.24	2' x 4'
T2262.33	3' x 3'
T2262.44	4' x 4'
T2262.46	4' x 6'
T2262.48	4' x 8'


Mini Plate Straight 2.5mm (Without Gap)




TT Gr.5	Size
T2263.02	2 Hole
T2263.03	3 Hole
T2263.04	4 Hole
T2263.05	5 Hole
T2263.06	6 Hole
T2263.08	8 Hole
T2263.10	10 Hole
T2263.12	12 Hole
T2263.16	16 Hole
T2263.20	20 Hole
T2263.30	30 Hole
T2263.50	50 Hole

Mini Plate Straight 2.5mm (With Gap)


TT Gr.5	Size
T2264.02	2 Hole
T2264.03	3 Hole
T2264.04	4 Hole
T2264.05	5 Hole
T2264.06	6 Hole
T2264.08	8 Hole



Mini 'L' Shape 90° Plate 2.5mm (Without Gap)





	TT Gr.5	Size
Right	T2265R.22	2+2 Hole
Left	T2265L.22	2+2 Hole

Mini 'L' Shape 90° Plate Right 2.5mm

TT Gr.5	Size
T2266R.22	2+2 Hole
T2266R.23	2+3 Hole
T2266R.32	3+2 Hole
T2266R.33	3+3 Hole

Mini 'L' Shape 90° Plate Left 2.5mm

TT Gr.5	Size
T2266L.22	2+2 Hole
T2266L.23	2+3 Hole
T2266L.32	3+2 Hole
T2266L.33	3+3 Hole



Left


Right

Mini 'L' Shape 100° Plate 2.5mm

	TT Gr.5	Size
Right	T2267R.32	3+2 Hole
Left	T2267L.32	3+2 Hole

Mini 'L' Shape 120° Plate 2.5mm(Without Gap)


	TT Gr.5	Size
Right	T2268R.22	2+2 Hole
Left	T2268L.22	2+2 Hole



Left

Right


Mini 'L' Shape 120° Plate Right 2.5mm




TT Gr.5	Size
T2269R.22	2+2 Hole
T2269R.23	2+3 Hole
T2269R.32	3+2 Hole
T2269R.33	3+3 Hole

Mini 'L' Shape 120° Plate Left 2.5mm

TT Gr.5	Size
T2269L.22	2+2 Hole
T2269L.23	2+3 Hole
T2269L.32	3+2 Hole
T2269L.33	3+3 Hole




Mini 'T' Shape Plate 2.5mm




TT Gr.5	Size
T2270.22	2+2 Hole
T2270.32	3+2 Hole
T2270.33	3+3 Hole

Mini 'Y' Shape Plate 2.5mm (Without Gap)

TT Gr.5	Size
T2271.22	2+2 Hole
T2271.23	2+3 Hole



Mini 'Y' Shape Plate 2.5mm (With Gap)



TT Gr.5	Size
T2272.22	2+2 Hole
T2272.23	2+3 Hole

Mini Double 'Y' Shape Plate 2.5mm



TT Gr.5

T2273.00


Mini 'I' Shape Plate 2.5mm

TT Gr.5

T2274.00



Mini 'H' Shape Plate 2.5mm



TT Gr.5

T2275.00


Mini 'X' Shape Plate 2.5mm

TT Gr.5

T2276.00



Mini Angular Plate 150° Right 2.5mm




TT Gr.5	Size
T2277R.04	4 Hole
T2277R.06	6 Hole

Mini Angular Plate 150° Left 2.5mm

TT Gr.5	Size
T2277L.04	4 Hole
T2277L.06	6 Hole




Mini Curved Plate Right 2.5mm




TT Gr.5	Size
T2278R.04	4 Hole
T2278R.06	6 Hole

Mini Curved Plate Left 2.5mm

TT Gr.5	Size
T2278L.04	4 Hole
T2278L.06	6 Hole





Mini 'Z' Shape Plate Type-1 2.5mm

TT Gr.5	Size
T2279R.00	Right
T2279L.00	Left

Mini 'Z' Shape Plate Type-2 Right 2.5mm



TT Gr.5	Size
T2280R.00	Right
T2280L.00	Left

**Mini Orbita Plate 2.5mm (With Gap)**

TT Gr.5	Size
T2281.04	4 Hole
T2281.06	6 Hole
T2281.08	8 Hole

Mini Orbita Plate 2.5mm (Without Gap)

TT Gr.5	Size
T2282.04	4 Hole
T2282.06	6 Hole
T2282.08	8 Hole
T2282.10	10 Hole
T2282.12	12 Hole
T2282.14	14 Hole
T2282.16	16 Hole
T2282.18	18 Hole
T2282.20	20 Hole

**Mini 3D Plate 2.5mm**

TT Gr.5	Size
T2283.22	2+2 Hole
T2283.23	2+3 Hole
T2283.24	2+4 Hole


Mandibular Reconstruction Plate Straight




TT Gr.5	Size
T2284.08	8 Hole
T2284.10	10 Hole
T2284.12	12 Hole
T2284.16	16 Hole
T2284.20	20 Hole
T2284.24	24 Hole

Mandibular Reconstruction Plate Double Angled 4 Hole

TT Gr.5	Size
T2285.16	16 Hole
T2285.20	20 Hole
T2285.24	24 Hole
T2285.28	28 Hole




Mandibular Reconstruction Plate Angled 4 Hole Right




TT Gr.5	Size
T2286R.12	12 Hole
T2286R.16	16 Hole
T2286R.20	20 Hole
T2286R.26	26 Hole

Mandibular Reconstruction Plate Angled 4 Hole Left

TT Gr.5	Size
T2286L.12	12 Hole
T2286L.16	16 Hole
T2286L.20	20 Hole
T2286L.26	26 Hole




Mandibular Reconstruction Plate Angled 6 Hole Right




TT Gr.5	Size
T2287R.16	16 Hole
T2287R.20	20 Hole

Mandibular Reconstruction Plate Angled 6 Hole Left


TT Gr.5	Size
T2287L.16	16 Hole
T2287L.20	20 Hole



Right




Left




TMJ Condyle Plate

	TT Gr.5	Size
Right	T2288R.03	3 Hole
Left	T2288L.03	3 Hole




TMJ Condyle With Reconstruction Plate 4 Hole Right




TT Gr.5	Size
T2289R.12	12 Hole
T2289R.16	16 Hole
T2289R.20	20 Hole
T2289R.24	24 Hole

TMJ Condyle With Reconstruction Plate 4 Hole left

TT Gr.5	Size
T2289L.12	12 Hole
T2289L.16	16 Hole
T2289L.20	20 Hole
T2289L.24	24 Hole




TMJ Condyle With Reconstruction Plate 6 Hole Right



TT Gr.5	Size
T2293R.16	16 Hole
T2293R.20	20 Hole

TMJ Condyle With Reconstruction Plate 6 Hole left


TT Gr.5	Size
T2293L.16	16 Hole
T2293L.20	20 Hole




Mandibular Reconstruction Plate U-Shape




TT Gr.5	Size
T2290.10	10 Hole
T2290.12	12 Hole
T2290.14	14 Hole
T2290.16	16 Hole

1.5mm Single Slotted Mini Screw


TT Gr.5	Size
T1080.004	4mm
T1080.005	5mm
T1080.006	6mm
T1080.007	7mm
T1080.008	8mm
T1080.009	9mm
T1080.010	10mm
T1080.011	11mm
T1080.012	12mm

1.5mm Cross Slotted Mini Screw

TT Gr.5	Size
T1081.004	4mm
T1081.005	5mm
T1081.006	6mm
T1081.007	7mm
T1081.008	8mm
T1081.009	9mm
T1081.010	10mm
T1081.011	11mm
T1081.012	12mm


1.7mm Single Slotted Emergency Screw

TT Gr.5	Size
T1082.006	6mm
T1082.007	7mm
T1082.008	8mm
T1082.009	9mm
T1082.010	10mm
T1082.011	11mm
T1082.012	12mm

1.7mm Cross Slotted Emergency Screw

TT Gr.5	Size
T1083.006	6mm
T1083.007	7mm
T1083.008	8mm
T1083.009	9mm
T1083.010	10mm
T1083.011	11mm
T1083.012	12mm


2.0mm Single Slotted Mini Screw




TT Gr.5	Size
T1084.004	4mm
T1084.005	5mm
T1084.006	6mm
T1084.007	7mm
T1084.008	8mm
T1084.009	9mm
T1084.010	10mm
T1084.011	11mm
T1084.012	12mm
T1084.013	13mm
T1084.014	14mm

2.0mm Cross Slotted Mini Screw

TT Gr.5	Size
T1085.004	4mm
T1085.005	5mm
T1085.006	6mm
T1085.007	7mm
T1085.008	8mm
T1085.009	9mm
T1085.010	10mm
T1085.011	11mm
T1085.012	12mm
T1085.013	13mm
T1085.014	14mm




2.3mm Single Slotted Emergency Screw




TT Gr.5	Size
T1086.004	4mm
T1086.005	5mm
T1086.006	6mm
T1086.007	7mm
T1086.008	8mm
T1086.009	9mm
T1086.010	10mm
T1086.011	11mm
T1086.012	12mm
T1086.013	13mm
T1086.014	14mm

2.3mm Cross Slotted Emergency Screw

TT Gr.5	Size
T1087.004	4mm
T1087.005	5mm
T1087.006	6mm
T1087.007	7mm
T1087.008	8mm
T1087.009	9mm
T1087.010	10mm
T1087.011	11mm
T1087.012	12mm
T1087.013	13mm
T1087.014	14mm



2.0mm Intermaxillary Fixation Screw




TT Gr.5	Size
T1088.004	4mm
T1088.006	6mm
T1088.008	8mm
T1088.010	10mm
T1088.012	12mm
T1088.014	14mm
T1088.016	16mm

2.0mm Leg Screw




TT Gr.5	Size
T1089.008	8mm
T1089.010	10mm
T1089.012	12mm
T1089.014	14mm
T1089.016	16mm
T1089.018	18mm
T1089.020	20mm

2.5mm Single Slotted Mini Screw




TT Gr.5	Size
T1090.004	4mm
T1090.005	5mm
T1090.006	6mm
T1090.007	7mm
T1090.008	8mm
T1090.009	9mm
T1090.010	10mm
T1090.011	11mm
T1090.012	12mm
T1090.013	13mm
T1090.014	14mm
T1090.015	15mm
T1090.016	16mm

2.5mm Cross Slotted Mini Screw




TT Gr.5	Size
T1091.004	4mm
T1091.005	5mm
T1091.006	6mm
T1091.007	7mm
T1091.008	8mm
T1091.009	9mm
T1091.010	10mm
T1091.011	11mm
T1091.012	12mm
T1091.013	13mm
T1091.014	14mm
T1091.015	15mm
T1091.016	16mm

2.7mm Single Slotted Emergency Screw




TT Gr.5	Size
T1092.004	4mm
T1092.005	5mm
T1092.006	6mm
T1092.007	7mm
T1092.008	8mm
T1092.009	9mm
T1092.010	10mm
T1092.011	11mm
T1092.012	12mm
T1092.013	13mm
T1092.014	14mm
T1092.015	15mm
T1092.016	16mm

2.7mm Cross Slotted Emergency Screw



TT Gr.5	Size
T1093.004	4mm
T1093.005	5mm
T1093.006	6mm
T1093.007	7mm
T1093.008	8mm
T1093.009	9mm
T1093.010	10mm
T1093.011	11mm
T1093.012	12mm
T1093.013	13mm
T1093.014	14mm
T1093.015	15mm
T1093.016	16mm

2.5mm Intermaxillary Fixation Screw






TT Gr.5	Size
T1094.004	4mm
T1094.006	6mm
T1094.008	8mm
T1094.010	10mm
T1094.012	12mm
T1094.014	14mm
T1094.016	16mm

2.5mm Leg Screw



TT Gr.5	Size
T1095.008	8mm
T1095.010	10mm
T1095.012	12mm
T1095.014	14mm
T1095.016	16mm
T1095.018	18mm
T1095.020	20mm

Our Products Range



SAMAY Surgical Pvt. Ltd.



Manufactured by....

Samay Surgical Pvt. Ltd

Avadh Industrial Park - A,
Survey No. 896 & 897,
Plot No. 12 & 13, Khambha - 360311,
Dist. Rajkot, Gujarat, India.

Contact Us

Call:- +91 78781 52154

E-mail:- info@samaysurgical.com

Web:- www.samaysurgical.com

Get access to our
complete catalog



Scan to Explore